



## BFA-500 Body fat scale

| Small Equipments

### Details Info

It can measure the body height, weight, fat and blood pressure after putting a proper coin. printing, pouncing the measuring result

#### Functions:

This machine is controlled by microcomputer. It can measure the body height, weight, fat and blood pressure after putting a proper coin .The measuring result can be printed and spoken out .Different kinds languages can be selected according to your requirement. The original Omron blood pressure meter is high quality and difficult to damage .The thermal printer makes the printing fast and clear, and this kind of printer is heavy –duty and easy to operate and maintain. The 6" high-clear LCD and touch screen make the operation quite easy.

#### Technical Data:

Height measurement range	100cm-200cm
Minimum division	0.5cm
Weight measurement range	8-200kg
Minimum division	0.1kg
Self-height	240cm
Self-weight	45kg