



CIA-900

| Immuno Assay Analyzer

☰ Details Info

Feature of this instrument

For individual testing and batch testing, Automatic timing and remind
 Touch-screen
 Add new projects without the software upgrade
 Store more than 50,000 messages and variety of ways to search
 Built-in quality control chart is automatically generated function
 Automatic identification of the project and batch number of the test card
 Including calibration card
 Support of barcode scanners, to input the information easily
 Compatible with the hospital information system

The main parameters of the instrument:

Display	7 inch touch screen
Test Speed	120 tests / hour
Interface	USB/RS232 / network interface
Storage capacity	50000
Size	220*60*130mm